

Clean D

Dental Aerosol Suction System

Designed to safely capture chairside aerosols and helps protect your operator from airborne contaminants.



Low Noise Level 65dB



1 HEPA + 2 Carbon Filter



22,000 Pa Suction Power



Medical Grade UV-C Lamp



Specifications

Model	Clean D
Power Supply (V)	230 VAC / 50, 60Hz
Electric Current (A)	13.0A
Output Power (kW)	1.3kW
Air Flow Rate	3.00 m ³ /min
Size WxDxH	370 x 340 x 1200 mm (not including arm)
Weight	Approx. 25kg
Noise Decibel	Approx. 65dB
Warranty	1 Year Manufacturing Defects

